

Pre-eclampsia diagnosis in 2nd and 3rd trimester of pregnancy

Improved pre-eclampsia diagnosis with sFlt-1 and PIGF • Prognosis of adverse outcome in women with suspected pre-eclampsia • Sensitive and reliable measurement on Thermo Scientific B·R·A·H·M·S KRYPTOR





Biomarkers for pre-eclampsia management Improving the diagnostic tools for pre-eclampsia evaluation

Pre-eclampsia is a progressive, pregnancy-related disorder with severe complications for mother and child. A timely diagnosis is needed in order to prevent maternal and fetal morbidity or mortality. In the absence of a specific therapy other than delivery the main objective of frequent patient monitoring is to detect deterioration of a patient's condition and to counteract maternal and fetal risk.





10% of pregnant women show unspecific signs and symptoms of pre-eclampsia



Only **one fifth of them** are actually developing pre-eclampsia¹

Diagnostic standard for pre-eclampsia ...



The "gold standard" for pre-eclampsia diagnosis – assessment of blood pressure and proteinuria – offers only a **poor sensitivity and specificity** with regards to origin of disease and prediction of maternal and perinatal outcome.^{2,3}



Serum sFlt-1 and PIGF determination adds significant clinical benefit to standard procedures

Determination of the biomarkers sFIt-1 (soluble FMS-like Tyrosine Kinase) and PIGF (Placental Growth Factor) in maternal blood have shown to significantly improve risk stratification among women presenting for pre-eclampsia evaluation. With the new high sensitivity assays Thermo Scientific[™] B·R·A·H·M·S[™] sFIt-1 KRYPTOR[™] and Thermo Scientific B·R·A·H·M·S PIGF plus KRYPTOR it is now possible to detect serum levels of both biomarkers reliably throughout pregnancy and thus improve pre-eclampsia management.

Clearly differentiating between pre-eclampsia

... significantly improved by sFlt-1 and PIGF serum measurement



Measuring sFIt-1 and PIGF starting in mid pregnancy in women with suspected pre-eclampsia significantly improves the current evaluation of patients – for a cost-effective patient management and improved care.

Pre-eclampsia diagnosis and prognosis of adverse outcome

The added value of sFlt-1 and PIGF

Improved diagnosis of pre-eclampsia with sFlt-1/PIGF ratio

Recent studies have proved the additional benefit of the sFIt-1/PIGF ratio in diagnosing pre-eclampsia:

- Measurement of sFIt-1 and PIGF levels in maternal serum, starting in mid pregnancy, can **confirm pre-eclampsia diagnosis**, with the sFIt-1/PIGF ratio having a superior diagnostic ability compared to either of the biomarkers alone.^{7,8}
- The addition of sFlt-1/PIGF ratio to Doppler ultrasound measurement improves the sensitivity and specificity in diagnosing pre-eclampsia compared to the Doppler measurement alone.⁸
- In women presenting with hypertension, the sFlt-1/ PIGF ratio is able to distinguish between those who will develop pre-eclampsia and those with chronic or gestational hypertension. Women with pre-eclampsia have a significantly higher sFlt-1/PIGF ratio than women with other hypertensive disorders or controls.^{4,5}

Therefore, the sFIt-1/PIGF ratio is a valuable tool for confirming or excluding the diagnosis of pre-eclampsia and offers a high clinical value for clinical management, counseling and risk anticipation.



PIGF and sFIt-1 were measured on KRYPTOR in parallel on samples from pregnant women with normal pregnancy outcome and patients with pre-eclampsia. At a cut-off of 85 for the sFIt-1/PIGF ratio, the sensitivity was calculated at 95% and the specificity at 84% for diagnosing pre-eclampsia.

The higher the sensitivity of a test the more women with pre-eclampsia are identified and can be advised for closer monitoring.

Figure 1 Improved pre-eclampsia diagnosis with sFlt-1/PIGF ratio ⁹



Prognosis of adverse outcome with sFIt-1/PIGF ratio

Rana et al. showed that **women with any subsequent adverse outcome** in addition to hypertension had a significantly higher sFIt-1/PIGF ratio than those women without, especially when presenting before week 34 (Figure 2).⁵



Women who needed to be delivered within the next

2 weeks after presentation had a significantly higher sFIt-1/ PIGF ratio than women who could continue with their pregnancy (Figure 3).⁵



Figure 2 Prediction of adverse outcome with sFlt-1/PlGF ratio in women presenting < 34 weeks' gestation⁵ **Figure 3** Prediction of duration of pregnancy with sFlt-1/PlGF ratio in women presenting < 34 weeks' gestation ⁵

The sFlt-1/PIGF ratio is also a potent predictor for subsequent maternal and fetal adverse outcome in women already diagnosed with pre-eclampsia and can support clinical decisions to avoid severe pregnancy complications.



The role of angiogenic factors Biomarker levels correlate with severity of disease

sFlt-1 and PIGF are counterparts

Although the cause of pre-eclampsia remains unclear, it is likely that the syndrome may be initiated by an imbalance of placental factors that induce endothelial dysfunction.



Figure 4 sFlt-1 acts as potent antagonist of PIGF and VEGF by adhering to the receptor-binding domains, thus preventing interaction with endothelial receptors and inducing endothelial dysfunction







Measuring maternal serum concentrations of sFlt-1 and PIGF can differentiate healthy women from women with pre-eclampsia.^{7,12} Changes in sFlt-1 and PIGF levels also appear to reflect the severity of the disease: early-onset pre-eclampsia is associated with greater changes in PIGF and sFlt-1 compared to late-onset pre-eclampsia.¹³

Pre-eclampsia management throughout pregnancy

Improving the outcomes for mother and child

PIGF and PAPP-A: First trimester screening for timely intervention

Combined screening for pre-eclampsia in weeks 10-13 can reliably identify women at risk for developing pre-eclampsia.

Combined first trimester screening includes

- serum PIGF and PAPP-A measurement,
- determination of mean arterial pressure (MAP), and
- Uterine Artery Pulsatility Index (UAPI)

resulting in a detection rate of ${>}90\%$ for a fixed false positive rate of 5%.14

An early identification of high-risk women allows for preventive measures and intensified monitoring. Administering low-dose aspirin (<150 mg/day) to high-risk women before 16 weeks of gestation can reduce the incidence of pre-eclampsia by 50%.¹⁵





Facts on pre-eclampsia

- Multisystem, life-threatening pregnancy-related disorder
- A main reason for maternal and fetal morbidity and mortality ^{16,17}
- Incidence: 2-8% of pregnancies
- Definition: New onset hypertension and proteinuria >20 weeks of gestation in previously normotensive women¹⁸
- HELLP (Hemolysis, Elevated Liver enzymes, Low Platelets): Severe pre-eclampsia variant occurring in ≈ 20% of symptomatic women; defined by additional affection of liver and coagulation system²⁰
- Eclampsia: Final stage of disease, associated with severe tonic-clonic seizures and coma as well as brain injury, cerebral edema and stroke²⁰



sFlt-1/PIGF ratio: Improved diagnosis and prognosis of adverse outcome

First symptoms of pre-eclampsia (hypertension, proteinuria) are observed after 20 weeks of gestation.²⁰

Diagnosis of pre-eclampsia is difficult, as pre-eclampsia can be confused with other diseases such as pregnancy-induced hypertension. By adding sFIt-1/PIGF ratio to the current diagnostic standard, the **diagnosis of pre-eclampsia** in a symptomatic woman can be confirmed or excluded.²

In women with diagnosed pre-eclampsia, the sFlt-1/PIGF ratio is a potent **predictor of subsequent maternal and fetal adverse outcome** and can be useful for further patient management.⁵

		Confirm or exlcude diagnosis of pre-eclampsia with sFit-1/PIGF ratio																				
		Prognosis of adverse outcome using sFIt-1/PIGF ratio																				
19	20)	21	22	23	24	25	26	27	28	29	30	31	32	33	34	4 35	36	37	' 38	39	40
		Early-onset (severe) pre-eclampsia											Intermediate- onset (medium) pre-eclampsia					noder- npsia				
		Severe and rapidly progressing form of pre-eclampsia with multiple complications and the need to deliver the baby preterm ²¹										Milder forms of pre-eclamp less complications occurrin 34 weeks of gestation ²¹					mpsia v ring aft	vith er				

Figure 6 First clinical symptoms of pre-eclampsia are observed >20 weeks of gestation. The gestational age at onset correlates with the severity of maternal and fetal consequences.²



Complete pre-eclampsia portfolio

From safe screening to improved diagnosis with B·R·A·H·M·S sFlt-1 and PIGF

Thermo Scientific B·R·A·H·M·S sFlt-1 KRYPTOR

Automated immunofluorescent assay for the quantitative determination of the concentration of sFlt-1 (soluble FMS-like Tyrosine Kinase 1, also known as VEGF receptor-1) in human serum.

- 75 determinations per kit
- 9 min incubation time
- Two-point calibration
- Monoparametric control kit, 3 levels
- Wide measuring range: 22-90000 pg/mL
- Excellent precision

With the lower and upper detection limits of 22 and 90000 pg/mL B·R·A·H·M·S sFlt-1 KRYPTOR provides the measuring range needed for a **reliable detection of clinical sFlt-1 values throughout pregnancy**.









Thermo Scientific B·R·A·H·M·S PIGF plus KRYPTOR

Automated immunofluorescent assay for the quantitative determination of the concentration of PIGF (Placental Growth Factor) in human serum. The assay is specific for the measurement of human free PIGF-1.

- 75 determinations per kit
- 29 min incubation time
- Single-point calibration
- Monoparametric control kit, 3 levels
- Wide measuring range: 3.6-7000 pg/mL
- Excellent precision

With a detection limit of 3.6 pg/mL and an upper limit of 7000 pg/mL B·R·A·H·M·S PIGF plus KRYPTOR provides the high sensitivity needed for **measuring PIGF levels in first trimester** as well as a wide measuring range to **reliably measure clinical values in second and third trimester**.

Exceptionally precise, fast and easy

Thermo Scientific B·R·A·H·M·S KRYPTOR compact PLUS

15 Years Reliable Results 15 Years Confident Decisions

- All KRYPTOR platforms FMF approved
- In routine use in labs worldwide since 1999
- Excellent precision and data stability



Thermo Scientific B·R·A·H·M·S Biomarkers Prenatal Screening Markers on KRYPTOR Systems

• BIRIAIHIMIS UE3 KRYPTOR	(under development)
 B·R·A·H·M·S sFlt-1 KRYPTOR* 	Art. no.: 845.075
 B·R·A·H·M·S PIGF plus KRYPTOR* 	Art. no.: 859.075
 B·R·A·H·M·S PIGF KRYPTOR* 	Art. no.: 844.075
• B·R·A·H·M·S PAPP-A KRYPTOR	Art. no.: 866.075
• B·R·A·H·M·S Inhibin A KRYPTOR	(under development)
 B·R·A·H·M·S hCG+β KRYPTOR 	Art. no.: 841.050
 B·R·A·H·M·S Free βhCG KRYPTOR 	Art. no.: 809.075
 B·R·A·H·M·S AFP KRYPTOR 	Art. no.: 816.075

• B·R·A·H·M·S Fast Screen pre I plus Software Art. no.: 105750

* Available on KRYPTOR compact PLUS

References

- 1. Milne et al. BMJ 2009; 339: b3129
- 2. Verlohren et al.Clin Sci 2012; 122(2): 43-52
- 3. Zhang et al. Obstet Gyneol 2001; 97: 261-7
- 4. Verlohren et al. Am J Obstet Gynecol 2012; 206, 58.e1-8
- 5. Rana et al. Circulation 2012; 125 (7): 911-9
- 6. Schnettler et al. BJOG 2013; 120: 1224-32
- 7. Verlohren et al. Am J Obstet Gynecol 2010; 202: 161.e1-11
- 8. Hagmann et al. Clin Chem 2012; 58(5), 837-45
- 9. B·R·A·H·M·S sFIt-1 KRYPTOR Instructions for use, 2013
- 10. Levine et al. N Engl J Med 2004; 350: 672-83
- 11. De Vivo et al. Acta Obstet et Gynecol 2008; 87: 837-42

- 12. Thadhani et al. J Clin Endo 2004; 89(2): 770-5
- 13. Romero et al. J Matern fetal Neonatal Med 2008; 21:9-23
- 14. Akolekar et al. Prenat Diagn 2011; 31: 66-74
- 15. Bujold et al. Obstet Gynecol 2010; 116: 402-14
- 16. WHO. World Health Report 2005: 63
- 17. Duley et al. Semin Perinatol 2009; 33(3): 130-7
- 18. Definition of the American College of Obstetrics and Gynecology
- 19. Sibai et al. Am J Obstet Gynecol 2009; 200: e481-7
 - 20. Powe et al. Circulation 2011;123: 2856-69
 - 21. Poon et al. Hypertension 2009; 53: 812-8

thermoscientific.com/brahms

© 2014 Thermo Fisher Scientific Inc. All rights reserved.

KRYPTOR is a registered trademark of CIS bio international, licensed for use by B·R·A·H·M·S, a part of Thermo Fisher Scientific. All other trademarks are the property of Thermo Fisher Scientific and its subsidiaries. All data regarding specifications, terms and pricing correspond to the existing knowledge at the time of the printing. We are not responsible for any errors, misprints or changes. Reprint, also in parts, solely with prior written consent of B·R·A·H·M·S GmbH.

Thermo Fisher Scientific products are distributed worldwide; not all intended uses and applications mentioned in this printing are registered in every country.

Clinical Diagnostics

Thermo Fisher Scientific B·R·A·H·M·S GmbH Neuendorfstr. 25 16761 Hennigsdorf Germany +49 (0)3302 883 0 +49 (0)3302 883 100 fax info.brahms@thermofisher.com www.thermoscientific.com/brahms www.thermoscientific.com/copeptin www.thermoscientific.com/proadrenomedullin www.thermoscientific.com/procalcitonin www.thermoscientific.com/kryptor



